

# Orbitor<sup>2</sup>

Pilling and snagging tester

**We specialise  
in instruments  
for evaluating  
the potential  
of woven or  
knitted fabrics  
to pill or snag.**

Whether you need to comply with a standard or a test method or to replicate a particular surface defect, we have the know-how to find the right solution for you.

Think of Orbitor<sup>2</sup> as a flexible framework, to which you can add different test boxes, drums or chambers, depending on the type of testing work you wish to carry out.

## Two models

To suit different volumes of testing, we offer two- or four-chamber instruments, with any combination of interchangeable test chambers.

Specimens are mounted on tubes and tumbled for an agreed period.



**Four-Position Orbitor<sup>2</sup>**

## Evaluation

After tumbling, the change in surface appearance is visually assessed under controlled viewing conditions.

For the evaluation of snagging, pinned bars or individual snagging points are added to the chambers.